

THE CARRIERE® SELF-LIGATING BRACKET SYSTEM

Shifting the way you think about orthodontics.



fast

natural

gentle



The Carriere Philosophy.

Fast.

Compared with conventional systems, the *Carriere* Self-Ligating Bracket System is faster, requires fewer appointments and is more comfortable for your patients. Quick and easy archwire changes and increased time between appointments afford a significant reduction in the number of patients seen per day, increasing practice profitability. Now, active treatment can be reduced by an average of four months.

Natural.

With the *Carriere* Self-Ligating Bracket, teeth are moved naturally with a low-force and low-friction passive system. This “freedom of fit” design stimulates cellular activity without occluding blood vessels in the periodontal ligament, allowing the teeth to actually move more quickly. These light forces are more biologically compatible than traditional, high force brackets and allow for more efficient tooth movement and improved patient comfort.

Gentle.

The *Carriere* Self-Ligating Bracket System produces significantly less friction than conventional or active self-ligating appliances. By greatly reducing the force and the pressure, low-force archwires can work to peak expression, thereby stimulating a more biologically compatible tooth movement. This unique low-force approach throughout the treatment means substantially greater patient comfort.

*Harradine NWT: Self-ligating brackets and treatment efficiency, 2001.

The Carriere Self-Ligating Bracket System

Takes clinical performance and patient outcomes to a new level.

The Carriere Self-Ligating Bracket is one of the fastest growing self-ligating brackets on the market today. It is specifically designed to meet your goal of low friction and low force by providing a “freedom of fit” in the bracket and archwire interface:

- In the body of the bracket, the mesial and distal edges of the slot have been carefully rounded for free sliding.
- The bracket archwire interface has a passive four-wall design that, with low super-elastic archwires, provides a free but controlled, synergistic action.

The “passive” orthodontic treatment philosophy can be “activated” when there is interest in using more force by using larger rectangular archwires.

Dual-Lock Fasteners in the Bracket Face

Cap opens to the incisal/occlusal edge allowing the cap to always stay in the closed position (e.g. during the mastication movements).

Smooth Labial Surface of the Locking Mechanism

Contoured edges and smooth surfaces provide enhanced patient comfort.

Simple Locking Mechanism

Opens with the Wire Director and Opener Tool and closes securely with finger.

Hooks Available on 3s, 4s and 5s

For maximum versatility.

Torque-In-Base/Compound-Contour

For ultimate precision, control and accurate fit.

Beveled Edges

Mesial and distal edges of the slot have been carefully rounded to reduce friction and enhance sliding mechanics.

Tie Wings

For chain elastics or conventional ligation (if necessary).



The Carriere Self-Ligating Bracket System features:

- A simple locking mechanism located in the bracket face, which can be easily opened with an explorer and closed with the touch of finger – offering quick, easy archwire changes.
- A familiar rhomboid shape that allows you to keep your preferred treatment techniques.
- Available in .018 and .022 McLaughlin, Bennett and Trevisi* and .018 and .022 Roth* prescriptions, available with or without hooks.

*Prescription does not imply endorsement by the doctors.




Fig. 01 | The Carriere with the LX Difference: Micro-etched surface for increased bond strength. Inset is 100x SEM.

Proprietary Improvements for a Stronger, Easier and Greater Performance.

The Carriere Self-Ligating Bracket with the LX difference will provide you with the following:

- Micro-etched surface for increased bond strength (Fig. 01).
- Lowest profile and less bulk in comparison to Damon 2 & 3™ and In-Ovation-R® (Fig. 02).
- Color-coded brackets for easier identification.
- Nickel free for increased biocompatibility and acceptance.
- Wire Director and Opener Tool – facilitates easier archwire placement and assists with opening and closing of the gate, REF 201-504W (Fig. 03).




Fig. 02 | SLB Profile Comparison

Accessories



Fig. 03 | Carriere Wire Director and Opener Tool 201-504W



Fig. 04 | Opti-Mim® Bite Corrector
Stainless Steel 430-006
Composite 430-007



Fig. 05 | Rotation Wedge
Silver 400-302
Clear 400-303

A Simplified Clinical Procedure.



Step 01 | Opening Locking Mechanism
Place tip of Wire Director and Opener (or instrument such as an explorer) in opening and apply light incisal/occlusal pressure.



Step 02 | Seating of the Archwire
Use the forked end of a Wire Director and Opener to help seat wire. Fully engage wire before closing the locking mechanism.



Step 03 | Closing Locking Mechanism
Close the locking mechanism by pressing the locking mechanism gingivally with your finger.

Carriere Bracket Sets

| | .018 | .022 | | .018 | .022 |
|----------------------|-----------|-----------|-----------------------------------|-----------|-----------|
| ROTH Rx* | | | MCLAUGHLIN, BENNETT & TREVISI Rx* | | |
| NO HOOKS | 703-297RQ | 703-397RQ | NO HOOKS | 703-297MQ | 703-397MQ |
| HOOKS ON 3S | 703-298RQ | 703-398RQ | HOOKS ON 3S | 703-298MQ | 703-398MQ |
| HOOKS ON 3S, 4S & 5S | 702-299RQ | 702-399RQ | HOOKS ON 3S, 4S & 5S | 703-299MQ | 703-399MQ |

Carriere McLaughlin, Bennett and Trevisi Rx*

| | COLOR | TOOTH # | TORQUE | ANGLE | .018 R | .018 L | .022 R | .022 L |
|-----------------------|--------|---------|--------|-------|----------------------|----------------------|----------------------|----------------------|
| MAXILLARY | | | | | | | | |
| CENTRALS | BLACK | U1 | +17° | +4° | 703-209 | 703-210 | 703-309 | 703-310 |
| LATERALS | PINK | U2 | +10° | +8° | 703-227 | 703-228 | 703-327 | 703-328 |
| CUSPIDS | GREEN | U3 | -7° | +8° | 703-253 703-253HK | 703-254 703-254HK | 703-353 703-353HK | 703-354 703-354HK |
| BICUSPIDS - 1ST & 2ND | PURPLE | U4&5 | -7° | 0° | 703-259 703-259HK | 703-260 703-260HK | 703-359 703-359HK | 703-360 703-360HK |
| MANDIBULAR | | | | | | | | |
| ANTERIORS | YELLOW | L1&2 | -6° | 0° | 703-237 | 703-238 | 703-337 | 703-338 |
| CUSPIDS | BLUE | L3 | -6° | +3° | 703-271 703-271HK | 703-272 703-272HK | 703-371 703-371HK | 703-372 703-372HK |
| BICUSPIDS - 1ST | GREY | L4 | -12° | 0° | 703-287 703-287HK | 703-288 703-288HK | 703-387 703-387HK | 703-388 703-388HK |
| BICUSPIDS - 2ND | WHITE | L5 | -17° | 0° | 703-289 703-289HK | 703-290 703-290HK | 703-389 703-389HK | 703-390 703-390HK |

Carriere Roth Rx*

| | COLOR | TOOTH # | TORQUE | ANGLE | .018 R | .018 L | .022 R | .022 L |
|-----------------------|--------|---------|--------|-------|----------------------|----------------------|----------------------|----------------------|
| MAXILLARY | | | | | | | | |
| CENTRALS | BLACK | U1 | +12° | +5° | 703-207 | 703-208 | 703-307 | 703-308 |
| LATERALS | PINK | U2 | +8° | +9° | 703-223 | 703-224 | 703-323 | 703-324 |
| CUSPIDS | GREEN | U3 | -2° | +11° | 703-251 703-251HK | 703-252 703-252HK | 703-351 703-351HK | 703-352 703-352HK |
| BICUSPIDS - 1ST & 2ND | PURPLE | U4&5 | -7° | 0° | 703-259 703-259HK | 703-260 703-260HK | 703-359 703-359HK | 703-360 703-360HK |
| MANDIBULAR | | | | | | | | |
| ANTERIORES | YELLOW | L1&2 | 0° | 0° | 703-235 | 703-236 | 703-335 | 703-336 |
| CUSPIDS | BLUE | L3 | -11° | +5° | 703-267 703-267HK | 703-268 703-268HK | 703-367 703-367HK | 703-368 703-368HK |
| BICUSPIDS - 1ST | WHITE | L4 | -17° | 0° | 703-289 703-289HK | 703-290 703-290HK | 703-389 703-389HK | 703-390 703-390HK |
| BICUSPIDS - 2ND | RED | L5 | -22° | 0° | 703-293 703-293HK | 703-294 703-294HK | 703-393 703-393HK | 703-394 703-394HK |

Carriere Self-Ligating Bracket System Recommended Archwire Sequences

| .018 | Pro Form™ Arch Form | | D-LX™ Form |
|--|---------------------|-----------|------------|
| | UPPER | LOWER | UNIVERSAL* |
| INITIAL | FIG. 01 | | FIG. 02 |
| .014 Dimpled <i>Bio-Kinetix</i> ® Thermally Activated <i>Nitanium</i> ® Archwire, <i>Pro Form</i> Arch Form | 100-852DM | 100-853DM | - |
| .014 <i>Bio-Kinetix Plus</i> Thermally Activated <i>Nitanium</i> Archwire, <i>D-LX</i> Arch Form | - | - | 103-101 |
| INTERMEDIATE | | | |
| .014 x .025 Dimpled <i>Bio-Kinetix</i> Thermally Activated <i>Nitanium</i> Archwire, <i>Pro Form</i> Arch Form | 100-872DM | 100-873DM | - |
| .014 x .025 <i>Bio-Kinetix Plus</i> Thermally Activated <i>Nitanium</i> Archwire, <i>D-LX</i> Arch Form | - | - | 103-104 |
| .016 x .025 Dimpled <i>Bio-Kinetix</i> Thermally Activated <i>Nitanium</i> Archwire, <i>Pro Form</i> Arch Form | 100-874DM | 100-875DM | - |
| .016 x .025 <i>Bio-Kinetix Plus</i> Thermally Activated <i>Nitanium</i> Archwire, <i>D-LX</i> Arch Form | - | - | 103-105 |
| FINISHING | | | |
| .016 x .022 <i>CNA</i> ® Beta III Archwire, <i>Pro Form</i> Arch Form | 100-944 | 100-954 | - |
| .016 X .025 <i>CNA</i> Beta III Archwire, <i>D-LX</i> Arch Form | - | - | 103-118 |
| .017 x .025 <i>CNA</i> Beta III Archwire, <i>Pro Form</i> Arch Form | 100-945 | 100-955 | - |
| .017 x .025 <i>CNA</i> Beta III Archwire, <i>D-LX</i> Arch Form | - | - | 103-119 |

Carriere Self-Ligating Bracket System Recommended Archwire Sequences (continued)




Fig. 01 | Pro Form Arch Form



Fig. 02 | D-LX Arch Form

| .022 | Pro Form Arch Form | | D-LX Form |
|---|--------------------|-----------|------------|
| | UPPER | LOWER | UNIVERSAL* |
| INITIAL | FIG. 01 | | FIG. 02 |
| .014 Dimpled <i>Bio-Kinetix</i> Thermally Activated Nitanium Archwire, <i>Pro Form Arch Form</i> | 100-852DM | 100-853DM | - |
| .014 <i>Bio-Kinetix Plus</i> Thermally Activated Nitanium Archwire, <i>D-LX Arch Form</i> | - | - | 103-101 |
| .016 Dimpled <i>Bio-Kinetix</i> Thermally Activated Nitanium Archwire, <i>Pro Form Arch Form</i> | 100-854DM | 100-855DM | - |
| .016 <i>Bio-Kinetix Plus</i> Thermally Activated Nitanium Archwire, <i>D-LX Arch Form</i> | - | - | 103-102 |
| OR | | | |
| .012 Dimpled Super Elastic Nitanium <i>Black-Ti</i> ® Archwire, <i>Pro Form Arch Form</i> | 101-901 | 101-902 | - |
| .014 Dimpled Super Elastic Nitanium <i>Black-Ti</i> Archwire, <i>Pro Form Arch Form</i> | 101-903 | 101-904 | - |
| INTERMEDIATE | | | |
| .014 x .025 Dimpled <i>Bio-Kinetix</i> Thermally Activated Nitanium Archwire, <i>Pro Form Arch Form</i> | 100-872DM | 100-873DM | - |
| .014 x .025 <i>Bio-Kinetix Plus</i> Thermally Activated Nitanium Archwire, <i>D-LX Arch Form</i> | - | - | 103-104 |
| .016 x .025 Dimpled <i>Bio-Kinetix</i> Thermally Activated Nitanium Archwire, <i>Pro Form Arch Form</i> | 100-874DM | 100-875DM | - |
| .016 x .025 <i>Bio-Kinetix Plus</i> Thermally Activated Nitanium Archwire, <i>D-LX Arch Form</i> | - | - | 103-105 |
| .018 x .025 Dimpled <i>Bio-Kinetix</i> Thermally Activated Nitanium Archwire, <i>Pro Form Arch Form</i> | 100-866DM | 100-867DM | - |
| .018 x .025 <i>Bio-Kinetix Plus</i> Thermally Activated Nitanium Archwire, <i>D-LX Arch Form</i> | - | - | 103-106 |
| FINISHING | | | |
| .017 x .025 CNA Beta III Archwire, <i>Pro Form Arch Form</i> | 100-945 | 100-955 | - |
| .017 x .025 CNA Beta III Archwire, <i>D-LX Arch Form</i> | - | - | 103-119 |
| .019 x .025 CNA Beta III Archwire, <i>Pro Form Arch Form</i> | 100-948 | 100-958 | - |
| .019 x .025 CNA Beta III Archwire, <i>D-LX Arch Form</i> | - | - | 103-120 |

*D-LX arch is offered exclusively in single packs. For additional information on our archwire selections and recommendations, please see our additional product brochure that can be downloaded at OrthoOrganizers.com or CarriereSystem.com, or contact your customer representative at 888.851.0533.

The Carriere Distalizer™ Appliance.



Also available for use with the *Carriere* Self-Ligating Bracket System. With its sleek, non-invasive design, the *Carriere Distalizer* Class II Correction Appliance provides greater comfort for your patients and shortens treatment time by up to 4 months. Delivering a natural, gentle and uniform force for distal molar movement with controlled rotation and tipping correction before brackets or other appliances are placed, the *Carriere Distalizer* Appliance easily turns your once complex Class II cases into simple Class I cases without extractions.

For more information on our products and educational offerings, please contact us:

In the U.S. 888.851.0533 | Outside the U.S. +(1) 760 448 8600 | Canada: CERUM 800.661.9567

To fax an order: 800.888.7244 | To email an order: USASales@OrthoOrganizers.com

CarriereSystem.com

© 2010 Ortho Organizers, Inc. All rights reserved. PN 999-205 Rev. 03/10. U.S. Patent No. 7,621,743 and foreign patents. Damon 2, Damon 3, and In-Ovation-R are the property of their respective owners. Ortho Organizers prescription does not imply endorsement by the doctor(s).

 **CARRIERE®**